

SR protein repeat Conjugated Antibody

Catalog No: #C57009

Package Size: #C57009-AF350 100ul #C57009-AF405 100ul #C57009-AF488 100ul #C57009-AF555 100ul #C57009-AF594 100ul
#C57009-AF647 100ul #C57009-AF680 100ul #C57009-AF750 100ul #C57009-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SR protein repeat Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	SR protein repeat Antibody detects endogenous levels of total SR protein repeat
Immunogen Description	A synthesized peptide derived from human SR protein repeat
Conjugates	AF Dye
Other Names	B52; SFRS6; Splicing factor; SRP55; SRSF6;
Accession No.	Uniprot:Q13247
Calculated MW	38-45kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.