

# C1s Conjugated Antibody

Catalog No: #C57012

Package Size: #C57012-AF350 100ul #C57012-AF405 100ul #C57012-AF488 100ul #C57012-AF555 100ul #C57012-AF594 100ul  
#C57012-AF647 100ul #C57012-AF680 100ul #C57012-AF750 100ul #C57012-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                     |
|-----------------------|-----------------------------------------------------|
| Product Name          | C1s Conjugated Antibody                             |
| Host Species          | Rabbit                                              |
| Clonality             | Monoclonal                                          |
| Isotype               | Rabbit IgG                                          |
| Purification          | Affinity-chromatography                             |
| Species Reactivity    | Human                                               |
| Specificity           | C1s Antibody detects endogenous levels of total C1s |
| Immunogen Description | A synthesized peptide derived from human C1s        |
| Conjugates            | AF Dye                                              |
| Other Names           | C1 esterase; C1S;                                   |
| Accession No.         | Uniprot:P09871                                      |
| Calculated MW         | 90kDa                                               |
| Storage               | Store at 4°C in dark for 6 months                   |

## Application Details

### Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.