

## ALS2CR1 Conjugated Antibody

Catalog No: #C57038

Package Size: #C57038-AF350 100ul #C57038-AF405 100ul #C57038-AF488 100ul #C57038-AF555 100ul #C57038-AF594 100ul  
#C57038-AF647 100ul #C57038-AF680 100ul #C57038-AF750 100ul #C57038-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                             |
|-----------------------|-------------------------------------------------------------|
| Product Name          | ALS2CR1 Conjugated Antibody                                 |
| Host Species          | Rabbit                                                      |
| Clonality             | Monoclonal                                                  |
| Isotype               | Rabbit IgG                                                  |
| Purification          | Affinity-chromatography                                     |
| Species Reactivity    | Human                                                       |
| Specificity           | ALS2CR1 Antibody detects endogenous levels of total ALS2CR1 |
| Immunogen Description | A synthesized peptide derived from human ALS2CR1            |
| Conjugates            | AF Dye                                                      |
| Other Names           | ALS2CR1; Nif3l1;                                            |
| Accession No.         | Uniprot:Q9GZT8                                              |
| Calculated MW         | 38kDa                                                       |
| Storage               | Store at 4°C in dark for 6 months                           |

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.