

IL-33 Conjugated Antibody

Catalog No: #C57041

Package Size: #C57041-AF350 100ul #C57041-AF405 100ul #C57041-AF488 100ul #C57041-AF555 100ul #C57041-AF594 100ul
#C57041-AF647 100ul #C57041-AF680 100ul #C57041-AF750 100ul #C57041-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-33 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Mouse Rat
Specificity	IL-33 Antibody detects endogenous levels of total IL-33
Immunogen Description	A synthesized peptide derived from human IL-33
Conjugates	AF Dye
Other Names	DVS27; IL 1F11 ;IL 33; IL1F11; IL33; Interleukin 1 family member 11; Interleukin 33; Interleukin 33 precursor; Interleukin33; NFEHEV; NFHEV;
Accession No.	Uniprot:Q8BVZ5
Calculated MW	33kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.