

HNF-4-alpha Conjugated Antibody

Catalog No: #C57044

Package Size: #C57044-AF350 100ul #C57044-AF405 100ul #C57044-AF488 100ul #C57044-AF555 100ul #C57044-AF594 100ul
#C57044-AF647 100ul #C57044-AF680 100ul #C57044-AF750 100ul #C57044-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HNF-4-alpha Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	HNF-4-alpha Antibody detects endogenous levels of total HNF-4-alpha
Immunogen Description	A synthesized peptide derived from human HNF-4-alpha
Conjugates	AF Dye
Other Names	FRTS4; Hepatic nuclear factor 4 alpha; Hepatocyte nuclear factor 4 alpha; Hepatocyte nuclear factor 4; HNF 4 alpha; HNF 4; HNF4; HNF4A; HNF4a7; HNF4a8; HNF4a9; Hnf4alpha; HNF4alpha10/11/12; MODY 1; MODY; MODY1; NR2A1; NR2A21; TCF; TCF14; Tcf4; Transcription factor 14; Transcription factor HNF4;
Accession No.	Uniprot:P41235
Calculated MW	53kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.