

Thrombomodulin / CD141 Conjugated Antibody

Catalog No: #C57050

Package Size: #C57050-AF350 100ul #C57050-AF405 100ul #C57050-AF488 100ul #C57050-AF555 100ul #C57050-AF594 100ul
#C57050-AF647 100ul #C57050-AF680 100ul #C57050-AF750 100ul #C57050-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thrombomodulin / CD141 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Mouse Rat
Specificity	Thrombomodulin / CD141 Antibody detects endogenous levels of total Thrombomodulin / CD141
Immunogen Description	A synthesized peptide derived from human Thrombomodulin / CD141
Conjugates	AF Dye
Other Names	AHUS6; BDCA3; CD141; CD141 antigen; Fetomodulin; Thbd; THPH12; THRM; Thrombomodulin;
Accession No.	Uniprot:P15306
Calculated MW	105,70kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.