

Importin 9 / RANBP9 Conjugated Antibody

Catalog No: #C57071

Package Size: #C57071-AF350 100ul #C57071-AF405 100ul #C57071-AF488 100ul #C57071-AF555 100ul #C57071-AF594 100ul
#C57071-AF647 100ul #C57071-AF680 100ul #C57071-AF750 100ul #C57071-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Importin 9 / RANBP9 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	Importin 9 / RANBP9 Antibody detects endogenous levels of total Importin 9 / RANBP9
Immunogen Description	A synthesized peptide derived from human Importin 9 / RANBP9
Conjugates	AF Dye
Other Names	Imp9; Imp9a; Imp9b; Ipo9; RanBP9;
Accession No.	Uniprot:Q96P70
Calculated MW	140kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.