

TPMT Conjugated Antibody

Catalog No: #C57082

Package Size: #C57082-AF350 100ul #C57082-AF405 100ul #C57082-AF488 100ul #C57082-AF555 100ul #C57082-AF594 100ul
#C57082-AF647 100ul #C57082-AF680 100ul #C57082-AF750 100ul #C57082-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TPMT Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	TPMT Antibody detects endogenous levels of total TPMT
Immunogen Description	A synthesized peptide derived from human TPMT
Conjugates	AF Dye
Other Names	Thiopurine methyltransferase; TPMT;
Accession No.	Uniprot:P51580
Calculated MW	30kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.