

CD32A + CD32B + CD32C Conjugated Antibody

Catalog No: #C57085

Package Size: #C57085-AF350 100ul #C57085-AF405 100ul #C57085-AF488 100ul #C57085-AF555 100ul #C57085-AF594 100ul
 #C57085-AF647 100ul #C57085-AF680 100ul #C57085-AF750 100ul #C57085-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD32A + CD32B + CD32C Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	CD32A + CD32B + CD32C Antibody detects endogenous levels of total CD32A + CD32B + CD32C
Immunogen Description	A synthesized peptide derived from human CD32A + CD32B + CD32C
Conjugates	AF Dye
Other Names	CD32; CD32A; CD32B; CD32C; CDw32; FCG2; FcGR; FCGR2; FCGR2A; FCGR2A1; FCGR2B; FCGR2C; Fcgr3; Fcgr3a; FcgRII; FcRII a; FCRII; FcRII b; FcRII c; FcRIIC; IGFR2;
Accession No.	Uniprot:P12318/P31994/P31995
Calculated MW	37-43kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.