

# Pyrophosphatase 1 Conjugated Antibody

Catalog No: #C57094

Package Size: #C57094-AF350 100ul #C57094-AF405 100ul #C57094-AF488 100ul #C57094-AF555 100ul #C57094-AF594 100ul  
#C57094-AF647 100ul #C57094-AF680 100ul #C57094-AF750 100ul #C57094-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                 |
|-----------------------|---------------------------------------------------------------------------------|
| Product Name          | Pyrophosphatase 1 Conjugated Antibody                                           |
| Host Species          | Rabbit                                                                          |
| Clonality             | Monoclonal                                                                      |
| Isotype               | Rabbit IgG                                                                      |
| Purification          | Affinity-chromatography                                                         |
| Species Reactivity    | Human Mouse Rat                                                                 |
| Specificity           | Pyrophosphatase 1 Antibody detects endogenous levels of total Pyrophosphatase 1 |
| Immunogen Description | A synthesized peptide derived from human Pyrophosphatase 1                      |
| Conjugates            | AF Dye                                                                          |
| Other Names           | inorganic pyrophosphatase 1; IOPPP; PP1; PPA1; PPase; Pyp;                      |
| Accession No.         | Uniprot:Q15181                                                                  |
| Calculated MW         | 33kDa                                                                           |
| Storage               | Store at 4°C in dark for 6 months                                               |

## Application Details

### Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.