

# SFT Conjugated Antibody

Catalog No: #C57096

Package Size: #C57096-AF350 100ul #C57096-AF405 100ul #C57096-AF488 100ul #C57096-AF555 100ul #C57096-AF594 100ul  
#C57096-AF647 100ul #C57096-AF680 100ul #C57096-AF750 100ul #C57096-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                     |
|-----------------------|-----------------------------------------------------|
| Product Name          | SFT Conjugated Antibody                             |
| Host Species          | Rabbit                                              |
| Clonality             | Monoclonal                                          |
| Isotype               | Rabbit IgG                                          |
| Purification          | Affinity-chromatography                             |
| Species Reactivity    | Human Mouse Rat                                     |
| Specificity           | SFT Antibody detects endogenous levels of total SFT |
| Immunogen Description | A synthesized peptide derived from human SFT        |
| Conjugates            | AF Dye                                              |
| Other Names           | SFT; UBC4/5; Ubch5; UBCH5A; Ube2d1;                 |
| Accession No.         | Uniprot:P51668                                      |
| Calculated MW         | 17kDa                                               |
| Storage               | Store at 4°C in dark for 6 months                   |

## Application Details

### Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.