

BRAT1 Conjugated Antibody

Catalog No: #C57106

Package Size: #C57106-AF350 100ul #C57106-AF405 100ul #C57106-AF488 100ul #C57106-AF555 100ul #C57106-AF594 100ul
#C57106-AF647 100ul #C57106-AF680 100ul #C57106-AF750 100ul #C57106-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BRAT1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	BRAT1 Antibody detects endogenous levels of total BRAT1
Immunogen Description	A synthesized peptide derived from human BRAT1
Conjugates	AF Dye
Other Names	BAAT1; brat1; C7orf27; RMFSL;
Accession No.	Uniprot:Q6PJG6
Calculated MW	75kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.