

CLIC4 Conjugated Antibody

Catalog No: #C57118

Package Size: #C57118-AF350 100ul #C57118-AF405 100ul #C57118-AF488 100ul #C57118-AF555 100ul #C57118-AF594 100ul
#C57118-AF647 100ul #C57118-AF680 100ul #C57118-AF750 100ul #C57118-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CLIC4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	CLIC4 Antibody detects endogenous levels of total CLIC4
Immunogen Description	A synthesized peptide derived from human CLIC4
Conjugates	AF Dye
Other Names	Chloride intracellular channel 4; Clic4; CLIC4L; HUH1; MC3S5; mtCLIC; p64H1;
Accession No.	Uniprot:Q9Y696
Calculated MW	29kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.