

DNAJA4 Conjugated Antibody

Catalog No: #C57123

Package Size: #C57123-AF350 100ul #C57123-AF405 100ul #C57123-AF488 100ul #C57123-AF555 100ul #C57123-AF594 100ul
#C57123-AF647 100ul #C57123-AF680 100ul #C57123-AF750 100ul #C57123-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNAJA4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	DNAJA4 Antibody detects endogenous levels of total DNAJA4
Immunogen Description	A synthesized peptide derived from human DNAJA4
Conjugates	AF Dye
Other Names	DNAJA4; MST104; MSTP104; PRO1472;
Accession No.	Uniprot:Q8WW22
Calculated MW	45kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.