

# FAM50A Conjugated Antibody

Catalog No: #C57124



Package Size: #C57124-AF350 100ul #C57124-AF405 100ul #C57124-AF488 100ul #C57124-AF555 100ul #C57124-AF594 100ul  
#C57124-AF647 100ul #C57124-AF680 100ul #C57124-AF750 100ul #C57124-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                           |
|-----------------------|-----------------------------------------------------------|
| Product Name          | FAM50A Conjugated Antibody                                |
| Host Species          | Rabbit                                                    |
| Clonality             | Monoclonal                                                |
| Isotype               | Rabbit IgG                                                |
| Purification          | Affinity-chromatography                                   |
| Species Reactivity    | Human Mouse Rat                                           |
| Specificity           | FAM50A Antibody detects endogenous levels of total FAM50A |
| Immunogen Description | A synthesized peptide derived from human FAM50A           |
| Conjugates            | AF Dye                                                    |
| Other Names           | fam50a; HXC26; XAP5;                                      |
| Accession No.         | Uniprot:Q14320                                            |
| Calculated MW         | 40kDa                                                     |
| Storage               | Store at 4°C in dark for 6 months                         |

## Application Details

Suggested Dilution:

- AF350 conjugated: most applications: 1: 50 - 1: 250
- AF405 conjugated: most applications: 1: 50 - 1: 250
- AF488 conjugated: most applications: 1: 50 - 1: 250
- AF555 conjugated: most applications: 1: 50 - 1: 250
- AF594 conjugated: most applications: 1: 50 - 1: 250
- AF647 conjugated: most applications: 1: 50 - 1: 250
- AF680 conjugated: most applications: 1: 50 - 1: 250
- AF750 conjugated: most applications: 1: 50 - 1: 250
- Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.