

AKT1/2/3 Conjugated Antibody

Catalog No: #C57129

Package Size: #C57129-AF350 100ul #C57129-AF405 100ul #C57129-AF488 100ul #C57129-AF555 100ul #C57129-AF594 100ul
 #C57129-AF647 100ul #C57129-AF680 100ul #C57129-AF750 100ul #C57129-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	AKT1/2/3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	AKT1 + AKT2 + AKT3 Antibody detects endogenous levels of total AKT1 + AKT2 + AKT3
Conjugates	AF Dye
Other Names	PKB;RAC;PRKBA;PKB-ALPHA;RAC-ALPHA; AKT1 ; AKT2; AKT3; Akt1/2/3; pan Akt;
Accession No.	P31749/P31751/Q9Y243
Calculated MW	56kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

A synthesized peptide derived from human AKT1 + AKT2 + AKT3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.