

PSIP1 / LEDGF Conjugated Antibody

Catalog No: #C57139

Package Size: #C57139-AF350 100ul #C57139-AF405 100ul #C57139-AF488 100ul #C57139-AF555 100ul #C57139-AF594 100ul
 #C57139-AF647 100ul #C57139-AF680 100ul #C57139-AF750 100ul #C57139-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PSIP1 / LEDGF Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	PSIP1 / LEDGF Antibody detects endogenous levels of total PSIP1 / LEDGF
Conjugates	AF Dye
Other Names	DFS70; LEDGF; p52; p75; PAIP; Psp1; PSIP2;
Accession No.	O75475
Calculated MW	75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Transcriptional coactivator involved in neuroepithelial stem cell differentiation and neurogenesis. Involved in particular in lens epithelial cell gene regulation and stress responses.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.