

SEC62 Conjugated Antibody

Catalog No: #C57144

Package Size: #C57144-AF350 100ul #C57144-AF405 100ul #C57144-AF488 100ul #C57144-AF555 100ul #C57144-AF594 100ul
#C57144-AF647 100ul #C57144-AF680 100ul #C57144-AF750 100ul #C57144-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SEC62 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	SEC62 Antibody detects endogenous levels of total SEC62
Conjugates	AF Dye
Other Names	Mouse, Rat, Human
Accession No.	Q99442
Calculated MW	50kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Required for preprotein translocation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.