

ABL2 Conjugated Antibody

Catalog No: #C59235

Package Size: #C59235-AF350 100ul #C59235-AF405 100ul #C59235-AF488 100ul #C59235-AF555 100ul #C59235-AF594 100ul
#C59235-AF647 100ul #C59235-AF680 100ul #C59235-AF750 100ul #C59235-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ABL2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	23B9
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	ABL2 Antibody detects endogenous levels of total ABL2
Immunogen Description	A synthesized peptide derived from human ABL2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	ABL2
Other Names	ABL2; ABLL; Tyrosine kinase ARG; kinase Arg;
Accession No.	Uniprot:P42684
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C in dark for 6 months.

Application Details

WB 1:500~1:2000 FC 1:50

Background

ABL2 is a cytoplasmic tyrosine kinase which is closely related to but distinct from ABL1. The similarity of the proteins includes the tyrosine kinase domains and extends amino-terminal to include the SH2 and SH3 domains. ABL2 is expressed in both normal and tumor cells. The ABL2 gene product is expressed as two variants bearing different amino termini, both approximately 12-kb in length.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.