

CD32 Conjugated Antibody

Catalog No: #C59339

Package Size: #C59339-AF350 100ul #C59339-AF405 100ul #C59339-AF488 100ul #C59339-AF555 100ul #C59339-AF594 100ul
#C59339-AF647 100ul #C59339-AF680 100ul #C59339-AF750 100ul #C59339-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD32 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	25C8
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	CD32 Antibody detects endogenous levels of total CD32
Immunogen Description	A synthesized peptide derived from human CD32
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	CD32
Other Names	CD32; CD32A; CD32B; CD32C; FCG2; FcGR; FCGR2; FCGR2A; FCGR2B; FCGR2C; Fcgr3a; FCRII; IGFR2;
Accession No.	Uniprot:P12318
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C in dark for 6 months.

Application Details

WB 1:500~1:2000 IHC 1:50~1:200

Background

Binds to the Fc region of immunoglobulins gamma. Low affinity receptor. By binding to IgG it initiates cellular responses against pathogens and soluble antigens.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.