

P2Y6 Conjugated Antibody

Catalog No: #C59804

Package Size: #C59804-AF350 100ul #C59804-AF405 100ul #C59804-AF488 100ul #C59804-AF555 100ul #C59804-AF594 100ul
#C59804-AF647 100ul #C59804-AF680 100ul #C59804-AF750 100ul #C59804-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P2Y6 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	43E7
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human
Specificity	P2Y6 Antibody detects endogenous levels of total P2Y6
Immunogen Description	A synthesized peptide derived from human P2Y6
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Target Name	P2RY6
Other Names	P2ry6; P2Y6; P2Y6 receptor; PP2891;
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C in dark for 6 months.

Application Details

WB 1:1000~1:5000 IHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50

Background

Receptor for extracellular UDP > UTP > ATP. The activity of this receptor is mediated by G proteins which activate a phosphatidylinositol-calcium second messenger system.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.