

RELM beta antibody Biotin conjugated

Catalog No: #C62101B

Package Size: #C62101B 100ul

Description

Product Name	RELM beta antibody Biotin conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Species Reactivity	Human
Immunogen Type	Peptide
Immunogen Description	peptide derived from Human RELM beta..
Conjugates	Biotin
Target Name	RELM beta
Calculated MW	12kDa
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.

Background

Probable hormone

Note: This product is for in vitro research use only and is not intended for use in humans or animals.