

TCR alpha/beta Antibody FITC Conjugated

Catalog No: #C62164F

Package Size: #C62164F 100ul



Description

Product Name	TCR alpha/beta Antibody FITC Conjugated
Host Species	Mouse
Purification	The antibody was purified by affinity chromatography.
Species Reactivity	Human
Conjugates	FITC
Target Name	TRAC; TRBC1
Formulation	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Storage	Store at 4°C in dark for 6 months

Note: This product is for in vitro research use only and is not intended for use in humans or animals.