

phospho-FAK (Ser910) Polyclonal Antibody Bio Conjugated

Catalog No: #C90918Bio

Package Size: #C90918Bio 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	phospho-FAK (Ser910) Polyclonal Antibody Bio Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	QEI(p-S)PP
Conjugates	Biotin
Target Name	phospho-FAK (Ser910)
Other Names	phospho-FAK (Ser910)
Accession No.	Uniprot:Q05397Gene ID:5747
Calculated MW	116
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Published Papers

el at., Cancer-associated fibroblasts promote lymphatic metastasis in cholangiocarcinoma via the PDGF-BB/PDGFR-β mediated paracrine signaling network. In Aging Dis on 2024 Feb 1 by Jian Yan, Gang Xiao,et al..PMID:37307823, , (2024)
[PMID:37307823](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.