

# CD11c Polyclonal Antibody FITC Conjugated

Catalog No: #C91504FITC

Package Size: #C91504FITC 100ul

## Description

|                    |                                                                                 |
|--------------------|---------------------------------------------------------------------------------|
| Product Name       | CD11c Polyclonal Antibody FITC Conjugated                                       |
| Host Species       | Rabbit                                                                          |
| Clonality          | Polyclonal                                                                      |
| Purification       | affinity purified by Protein A                                                  |
| Applications       | Flow-Cyt ICC IF                                                                 |
| Species Reactivity | Hu Ms Rat                                                                       |
| Conjugates         | FITC                                                                            |
| Target Name        | CD11c                                                                           |
| Other Names        | CD11c                                                                           |
| Accession No.      | Uniprot:P20702Gene ID:3687                                                      |
| Calculated MW      | 126                                                                             |
| Concentration      | 1mg/ml                                                                          |
| Formulation        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

## Published Papers

Wang Yu-Han;Li A-Guo;Wang Hong-Yan;Tu Yong-Sheng et al., Repression of JAK2-STAT1 and PD-L1 by CEP-33779 ameliorates the LPS-induced decline in phagocytic activity of alveolar macrophages and mitigates lung injury in mice, , (2024)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.