

Salmonella Antibody HRP Conjugated

Catalog No: #C92912H

Package Size: #C92912H 100ul

Description

Product Name	Salmonella Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Species Reactivity	(predicted: Mouse, Salmonella)
Immunogen Type	Peptide
Immunogen Description	Whole cell protein of Salmonella
Conjugates	HRP
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.