

Fluorescein antibody HRP Conjugated

Catalog No: #C92917H

Package Size: #C92917H 100ul

Description

Product Name	Fluorescein antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF,ELISA
Species Reactivity	Human
Immunogen Description	KLH conjugated FITC
Conjugates	Unconjugated.
Target Name	FITC
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF 1:100-500;ELISA=1:2000-10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.