

VRK1 Antibody FITC Conjugated

Catalog No: #C92918F

Package Size: #C92918F 100ul

Description

Product Name	VRK1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB,IHC,IF
Species Reactivity	predicted: Human,Mouse,Rat,Pig,Sheep,Cow,Horse
Specificity	The antibody detects endogenous level of total VRK1 protein.
Immunogen Type	Protein
Immunogen Description	KLH conjugated synthetic peptide derived from human VRK1: 132-230/396
Conjugates	Unconjugated.
Target Name	VRK1
Calculated MW	45kD(Calculated)
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000;IHC 1:100-1:500;IF 1:100-1:500

Note: This product is for in vitro research use only and is not intended for use in humans or animals.