

IL12/Biotin Conjugated antibody

Catalog No: #C92919B

Package Size: #C92919B 100ul

Description

Product Name	IL12/Biotin Conjugated antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB,IHC,IFoΩ%oΩ%ELISA
Species Reactivity	Mouse, Rat
Specificity	The antibody detects endogenous level of total IL12 protein.
Immunogen Type	Protein
Immunogen Description	KLH conjugated synthetic peptide derived from mouse IL-12:51-150/215?
Conjugates	Biotin
Target Name	IL-12
Calculated MW	22kD
Concentration	1mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.