

HSV1 gD antibody FITC Conjugated

Catalog No: #C92925F

Package Size: #C92925F 100ul

Description

Product Name	HSV1 gD antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity purified by Protein A.
Applications	ICC IF
Immunogen Type	Peptide
Immunogen Description	peptide derived from human HSV1 gD, strain Angelotti
Conjugates	FITC
Target Name	HSV1 gD
Calculated MW	41kDa
Concentration	1mg/ml
Formulation	Liquid in 0.01M TBS (pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Store at +4°C for short term. Store at -20°C for long term. Avoid freeze/thaw cycle.

Application Details

ICC 1:50-200, IF 1:50-200

Background

Herpes simplex type 1 (HSV-1) belongs to a family that includes HSV-2, Epstein-Barr virus (EBV) and Varicella zoster (chicken pox) virus amongst others. HSV-1 and HSV-2 are extremely difficult to distinguish from each other. Members of this family have a characteristic virion structure. The double stranded DNA genome is contained within an icosahedral capsid embedded in a proteinaceous layer (tegument) and surrounded by a lipid envelope, derived from the nuclear membrane of the last host, which is decorated with virus-specific glycoproteins spikes. These viruses are capable of entering a latent phase where the host shows no visible sign of infection and levels of infectious agent become very low. During the latent phase the viral DNA is integrated into the genome of the host cell. Glycoprotein D (gD) has been implicated in binding to cellular receptors that facilitate virus penetration into cells. Herpes simplex virus type 1 (HSV-1) glycoprotein D (gD) is an essential component of the entry apparatus that is responsible for viral penetration and subsequent cell-cell spread.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.