

Hepatitis C Virus Core Antigen antibody Biotin conjugated

Catalog No: #C92938B

Package Size: #C92938B 100ul

Description

Product Name	Hepatitis C Virus Core Antigen antibody Biotin conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ELISA
Species Reactivity	Hepatitis C Virus
Immunogen Description	KLH conjugated synthetic peptide derived from human HCV-Core: 1-50/60
Calculated MW	19/21kDa
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ELISA 1:5000-10000

Background

HCV is classified into the genus Hepacivirus of the family Flaviviridae. Like all the members of the family, HCV is an enveloped, single-stranded, positive-sense RNA virus. Its genome (about 9600 nt) is flanked at both termini by conserved, highly structured non-translated regions (NTRs) and encodes a polyprotein precursor (about 3000 aa), which is proteolytically processed by host and viral proteases to produce the structural (core, E1, E2 and p7) and non-structural (NS2, NS3, NS4A, NS4B, NS5A and NS5B) proteins of the virus. Recently, an additional protein has been identified, whose function remains unknown. NS5A is a ~56 kDa pleiotropic protein with key roles in both viral RNA replication and modulation of the physiology of the host cell. It's exact role is not currently known (2008).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.