

Rat Cdk12/Cyclin-dependent kinase 12 ELISA Kit

Catalog No: #EK20000

Package Size: #EK20000 1x96T

Description

Product Name	Rat Cdk12/Cyclin-dependent kinase 12 ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cyclin-dependent kinase 12
Target Name	Cyclin-dependent kinase 12
Other Names	Cdc2-related kinase, arginine/serine-rich,CrkRS,Cell division cycle 2-related protein kinase 7,CDC2-related protein kinase 7,Cell division protein kinase 12,Protein kinase for splicing component,Crk7,Crkrs,Pksc
Accession No.	Q3MJK5

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:42pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cdk12

Note: This product is for in vitro research use only and is not intended for use in humans or animals.