

Mouse Cdh6/Cadherin-6 ELISA Kit

Catalog No: #EK20001

Package Size: #EK20001 1x96T



Description

Product Name	Mouse Cdh6/Cadherin-6 ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cadherin-6
Target Name	Cadherin-6
Other Names	Kidney cadherin,K-cadherin
Accession No.	P97326

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.087ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cdh6

Note: This product is for in vitro research use only and is not intended for use in humans or animals.