

Mouse Bin1/Myc box-dependent-interacting protein 1 ELISA Kit

Catalog No: #EK20024

Package Size: #EK20024 1x96T

Description

Product Name	Mouse Bin1/Myc box-dependent-interacting protein 1 ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Myc box-dependent-interacting protein 1
Target Name	Myc box-dependent-interacting protein 1
Other Names	Amphiphysin II,Amphiphysin-like protein,Bridging integrator 1,SH3 domain-containing protein 9,AmphI,Sh3p9
Accession No.	O08539

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.027ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Bin1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.