

Rat Chia/Acidic mammalian chitinase ELISA Kit

Catalog No: #EK20027

Package Size: #EK20027 1x96T



Description

Product Name	Rat Chia/Acidic mammalian chitinase ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Acidic mammalian chitinase
Target Name	Acidic mammalian chitinase
Other Names	AMCase
Accession No.	Q6RY07

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.075ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Chia

Note: This product is for in vitro research use only and is not intended for use in humans or animals.