

Mouse Avpr1b/Vasopressin V1b receptor ELISA Kit

Catalog No: #EK20028

Package Size: #EK20028 1x96T



Description

Product Name	Mouse Avpr1b/Vasopressin V1b receptor ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Vasopressin V1b receptor
Target Name	Vasopressin V1b receptor
Other Names	V1bR,AVPR V1b,AVPR V3,Antidiuretic hormone receptor 1b,Vasopressin V3 receptor,V1br
Accession No.	Q9WU02

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Avpr1b

Note: This product is for in vitro research use only and is not intended for use in humans or animals.