

Human ABCG1/ATP-binding cassette sub-family G member 1 ELISA Kit

Catalog No: #EK20030

Package Size: #EK20030 1x96T



Description

Product Name	Human ABCG1/ATP-binding cassette sub-family G member 1 ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human ATP-binding cassette sub-family G member 1
Target Name	ATP-binding cassette sub-family G member 1
Other Names	ATP-binding cassette transporter 8,White protein homolog,ABC8,WHT1
Accession No.	P45844

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.085ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:ABCG1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.