

Dog ANPEP/Aminopeptidase N ELISA Kit

Catalog No: #EK20032

Package Size: #EK20032 1x96T

Description

Product Name	Dog ANPEP/Aminopeptidase N ELISA Kit
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Aminopeptidase N
Target Name	Aminopeptidase N
Other Names	AP-N,cAPN,Alanyl aminopeptidase,Aminopeptidase M,AP-M,Microsomal aminopeptidase,CD13
Accession No.	P79143

Application Details

Assay Type: SandwichDetect Range:78-5000pg/ml Sensitivity:41.2pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:ANPEP

Note: This product is for in vitro research use only and is not intended for use in humans or animals.