

Human ANPEP/Aminopeptidase N ELISA Kit

Catalog No: #EK20033

Package Size: #EK20033 1x96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ANPEP/Aminopeptidase N ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Aminopeptidase N
Target Name	Aminopeptidase N
Other Names	AP-N, Alanyl aminopeptidase, Aminopeptidase M, AP-M, Microsomal aminopeptidase, Myeloid plasma membrane glycoprotein CD13, gp150, CD13, APN, CD13, PEPN
Accession No.	P15144

Application Details

Assay Type: Sandwich Detect Range: 7.8-500 ng/mL Sensitivity: 4.1 ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100 µl Assay Time: 3-5 h Detection wavelength: 450 nm gene_full_name: ANPEP

Note: This product is for in vitro research use only and is not intended for use in humans or animals.