

Mouse C8b/Complement component C8 beta chain ELISA Kit

Catalog No: #EK20035

Package Size: #EK20035 1x96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse C8b/Complement component C8 beta chain ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Complement component C8 beta chain
Target Name	Complement component C8 beta chain
Other Names	Complement component 8 subunit beta
Accession No.	Q8BH35

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.048ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: C8b

Note: This product is for in vitro research use only and is not intended for use in humans or animals.