

Human B3GAT1/Galactosylgalactosylxylosylprotein
3-beta-glucuronosyltransferase 1 ELISA Kit

Catalog No: #EK20037

Package Size: #EK20037 1x96T



Description

Product Name	Human B3GAT1/Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 1 ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 1
Target Name	Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 1
Other Names	Beta-1,3-glucuronyltransferase 1,Glucuronosyltransferase P,GlcAT-P,UDP-GlcUA:glycoprotein beta-1,3-glucuronyltransferase,GlcUAT-P,GLCATP
Accession No.	Q9P2W7

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:B3GAT1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.