

Human ADCY9/Adenylate cyclase type 9 ELISA Kit

Catalog No: #EK20039

Package Size: #EK20039 1x96T



Description

Product Name	Human ADCY9/Adenylate cyclase type 9 ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Adenylate cyclase type 9
Target Name	Adenylate cyclase type 9
Other Names	ATP pyrophosphate-lyase 9,Adenylate cyclase type IX,ACIX,Adenylyl cyclase 9,AC9,KIAA0520
Accession No.	O60503

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.172ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:ADCY9

Note: This product is for in vitro research use only and is not intended for use in humans or animals.