

Pig ATP2A3/Sarcoplasmic/endoplasmic reticulum calcium ATPase 3 ELISA Kit

Catalog No: #EK20040

Package Size: #EK20040 1x96T

Description

Product Name	Pig ATP2A3/Sarcoplasmic/endoplasmic reticulum calcium ATPase 3 ELISA Kit
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Sarcoplasmic/endoplasmic reticulum calcium ATPase 3
Target Name	Sarcoplasmic/endoplasmic reticulum calcium ATPase 3
Other Names	SERCA3, Calcium pump 3
Accession No.	O77696

Application Details

Assay Type: Sandwich Detect Range: 0.312-20 ng/mL Sensitivity: 0.173 ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100 µl Assay Time: 3-5 h Detection wavelength: 450 nm gene_full_name: ATP2A3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.