

Mouse Cd63/CD63 antigen ELISA Kit

Catalog No: #EK20042

Package Size: #EK20042 1x96T



Description

Product Name	Mouse Cd63/CD63 antigen ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse CD63 antigen
Target Name	CD63 antigen
Other Names	CD63
Accession No.	P41731

Application Details

Assay Type: SandwichDetect Range:31.2-2000pg/mL Sensitivity:15.6pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cd63

Note: This product is for in vitro research use only and is not intended for use in humans or animals.