

Rat Cd63/CD63 antigen ELISA Kit

Catalog No: #EK20043

Package Size: #EK20043 1x96T



Description

Product Name	Rat Cd63/CD63 antigen ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat CD63 antigen
Target Name	CD63 antigen
Other Names	Mast cell antigen AD1,CD63
Accession No.	P28648

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:39.7pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cd63

Note: This product is for in vitro research use only and is not intended for use in humans or animals.