

Mouse Aplnr/Apelin receptor ELISA Kit

Catalog No: #EK20044

Package Size: #EK20044 1x96T

Description

Product Name	Mouse Aplnr/Apelin receptor ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Apelin receptor
Target Name	Apelin receptor
Other Names	Angiotensin receptor-like 1,G-protein coupled receptor APJ,MSR,Agtr11,Apj
Accession No.	Q9WV08

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Aplnr

Note: This product is for in vitro research use only and is not intended for use in humans or animals.