

Mouse Bcam/Basal cell adhesion molecule ELISA Kit

Catalog No: #EK20045

Package Size: #EK20045 1x96T

Description

Product Name	Mouse Bcam/Basal cell adhesion molecule ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Basal cell adhesion molecule
Target Name	Basal cell adhesion molecule
Other Names	B-CAM cell surface glycoprotein,Lutheran antigen,CD239,Gplu,Lu
Accession No.	Q9R069

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.067ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Bcam

Note: This product is for in vitro research use only and is not intended for use in humans or animals.