

Rat Bcam/Basal cell adhesion molecule ELISA Kit

Catalog No: #EK20046

Package Size: #EK20046 1x96T



Description

Product Name	Rat Bcam/Basal cell adhesion molecule ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Basal cell adhesion molecule
Target Name	Basal cell adhesion molecule
Other Names	B-CAM cell surface glycoprotein,Lutheran antigen,CD239,Lu
Accession No.	Q9ESS6

Application Details

Assay Type: SandwichDetect Range:15.6-1000pg/mL Sensitivity:8.2pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Bcam

Note: This product is for in vitro research use only and is not intended for use in humans or animals.