

Mouse Cdkn1a/Cyclin-dependent kinase inhibitor 1 ELISA Kit

Catalog No: #EK20048

Package Size: #EK20048 1x96T

Description

Product Name	Mouse Cdkn1a/Cyclin-dependent kinase inhibitor 1 ELISA Kit
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cyclin-dependent kinase inhibitor 1
Target Name	Cyclin-dependent kinase inhibitor 1
Other Names	CDK-interacting protein 1,Melanoma differentiation-associated protein,p21,Cip1,Waf1
Accession No.	P39689

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:37pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cdkn1a

Note: This product is for in vitro research use only and is not intended for use in humans or animals.