

Rat Apoc4/Apolipoprotein C-IV ELISA Kit

Catalog No: #EK20050

Package Size: #EK20050 1x96T



Description

Product Name	Rat Apoc4/Apolipoprotein C-IV ELISA Kit
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Apolipoprotein C-IV
Target Name	Apolipoprotein C-IV
Other Names	Apo-CIV,Apolipoprotein C4,Apolipoprotein E-linked,ECL,Ecl
Accession No.	P55797

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:39pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Apoc4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.